

PLANILUX iris light 100, used for the accurate detection of minor contrasts, for the ...

??? ?????: 1001700

180,00 €
????? ??????



??? ??????

PLANILUX iris light 100, used for the accurate detection of minor contrasts, for the improved assessment of over-exposed X ray images with high light density, max. 70.000 CD/qm, continuous dimming, precision iris diaphragm of 5 - 100 mm in diameter, continuous collimation to the smallest image areas, in mint condition

-> ????? ??? ?????? ??? ?????? <-

