

**DR.MÜLLER extension arm for large animal circle section. This support arm is used to ...**

??? ?????: 1005310

**159,00 €**  
????? ??????



??? ??????

DR.MÜLLER extension arm for large animal circle section. This support arm is used to improve the manoeuvrability of the large animal circle section and to increase the distance to the base unit. Made of aluminium with 2 locking screws, dimensions: length 230mm, height 50mm, depth 34mm, new. This article is intended exclusively for use in veterinary medicine.

-> ????? ??? ?????? ??? ?????? <-

