

ERBE Erbotom T400, mono- and bipolar, 400 Watt, dimensions: HxWxD 17cmx48cmx34cm, ...

??? ?????: 1001306

1.390,00 €
????? ??????



??? ??????

ERBE Erbotom T400, mono- and bipolar, 400 Watt, dimensions: HxWxD 17cmx48cmx34cm, second-hand. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONO AND BIPOLAR ELECTRODES, 3 DIFFERENT INSERTS)

-> ????? ??? ?????? ??? ?????? <-

