

**SMITHS EQ-5000 Level 1 Equator,  
convective patient warmer. Adjustable  
temperature range: ...**

??? ?????: 1005227

**670,00 €**  
????? ??????



??? ??????

SMITHS EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: 36°C, 40°C, 44°C. Delivery incl. supply hose 2meter, used. Delivery includes 5-foot chassisWe recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: [vet@moeckundmoeck.de](mailto:vet@moeckundmoeck.de)

-> ????? ??? ?????? ??? ?????? <-

