

**GEUDER G-28252ECV Megatron Cool" III-PV" (version ext. sterilizable pressure transducer) ...**

??? ?????: 1004856

**2.900,00 €**  
????? ??????



??? ??????

GEUDER G-28252ECV Megatron Cool" III-PV" (version ext. sterilizable pressure transducer) ophthalmosurgical system, microprocessor-controlled unit with I/A system, ultrasound incl. Cool Flash" mode for MICS" (bimanualPako), vitrectomy (magnetic and pneumatic) and diathermy; incl. pressure transducer G-29850 With G-28540V combination foot switch FS4 Both units are tested by the manufacturer The devices offered here in the repository can be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. Please see our general terms and conditions for further information.

-> ????? ??? ??????? ??? ??????? <-

