

DRÄGER Perseus A500, mobile anaesthesia workstation for adults, children and neonates ...

Article number: 1005315

19.800,00 €
net price



Product description

DRÄGER Perseus A500, mobile anaesthesia workstation for adults, children and neonates Ventilation modes: MAN/SPON, PC-CMV, PC-BIPAP, VC-CMV, VC-SIMV, VC-CMV/AF, VC-SIMV/AF. Integrated emergency power supply, minimum 30 min, normally 150 min. To be operated with Dräger Vapor 2000 or 3000 (not included). Electronic TurboVent 2 fan. Pull-out writing tablet, integrated lighting, 2 drawers. Dimensions: (HxWxD: 148cmx115cmx79cm). Weight approx. 160kg, second-hand

-> [Click here to come to the site.](#) <-

