

NIHON KOHDEN BSM-3763 LifeScope, multi-parameter monitor for measuring ECG, SPO2, NIBP, ...

Article number: 1004907

2.980,00 €
net price



Product description

NIHON KOHDEN BSM-3763 LifeScope, multi-parameter monitor for measuring ECG, SPO2, NIBP, temperature, modular smart cable technology (The multi-connector socket automatically recognizes the parameter type: BIS, IBP, HZV, CO2). Extra large color touch display, network compatible, used

-> [Click here to come to the site.](#) <-

