

**SET with reusable tubing, NEW,  
consisting of: DRAEGER patient tube  
silicone 110 cm, 22 ...**

Numéro d'article: 1004628

**212,00 €**  
prix net



**Product description**

SET with reusable tubing, NEW, consisting of: DRAEGER patient tube silicone 110 cm, 22 mm, Y-piece angled with LUER Lock connection, DRAEGER breathing bag 2.0 liter made of chloroprene, DRAEGER iso-spout, DRAEGER patient tube made of silicone, 50 cm, 22 mm

-> [Cliquez ici pour entrer sur le site.](#) <-

